

# *Aneroid Sphygmomanometer*

## *BK2001*



### **Specification:**

<i>Product Name</i>	<i>Aneroid Sphygmomanometer</i>
<i>Model Number</i>	<i>BK2001</i>
<i>Measure scope</i>	<i>0-300mmHg</i>
<i>Accuracy</i>	<i>+/-3mmHg</i>
<i>Sub-division</i>	<i>2mmHg</i>
<i>Aneroid Gauge</i>	<i>110g/120g/140g</i>
<i>Cuff</i>	<i>023 Nylon/cotton cuff with/without D-ring 48X14cm/51X14cm</i>
<i>Bladder</i>	<i>PVC/Latex bladder with two tube, size 22X12cm</i>
<i>Carrying Bag</i>	<i>069 leather bag / 068 nylon bag</i>
<i>Stethoscope</i>	<i>Single head stethoscope</i>
<i>Chest piece</i>	<i>Aluminum chest piece</i>
<i>Y-tube</i>	<i>60cm long Y-tube</i>
<i>Color</i>	<i>black,blue,greys,red,green etc</i>