

desktop low-speed centrifuge 80-2



Features:

Desktop electric centrifuge 80-2 is a conventional laboratory centrifuge, widely used in biological, chemical, pharmaceutical and other research and education and production sectors. Low-speed desktop centrifuge maximum speed 4000R / MIN, suitable for radioimmunoassay and separation of cells and large particles, there are two kinds of timing and non-timing.

Technical Parameters:

Working power: 220V 50HZ

Maximum speed: 4000r / min (step speed)

Maximum relative centrifugal force: $1435 \times g$

Working capacity: $12 \times 20\text{ml}$

With time: 0-60 minutes

Model	Max. speed	Capacity of the rotor	Max RCF	Timer	Power
80-2	4000rpm	20ml×12	1795xg	0-60min or continually	AC220V 50HZ