

# DRY OVEN

## DHG-9023A

- Microprocessor control.
- With double glass window.
- Digital display setting and actual temperature.
- Digital display running time.
- PID self-tuning, over temperature alarm.
- Automatic alarm after operation.
- Stainless-steel interior, chamber corners are semicircle.
- It is applicable to drying, baking, sterilization, dewaxing, solidifying for industrial and mining enterprises, pharmaceutical and health care, and scientific research institutes.



MODEL	DHG-9023A
Power	320W
Working size(mm)	300x300x275
Temperature Range	50~300°C
Temperature Resolution	0.1°C
Temperature Fluctuation	±1°C
Time	0-999min
Voltage	110V 60Hz / 220v 50Hz
Dimensions(mm)	570x530x640
Weight	28kg