

HOT OVEN

DHG-9041A

- Microprocessor control.
- With double glass window.
- Digital display setting and actual temperature.
- Digital display running time.
- PID self-tuning, over temperature alarm.
- Automatic alarm after operation.
- Stainless-steel interior, chamber corners are semicircle.
- It is applicable to drying, baking, sterilization, dewaxing, solidifying for industrial and mining enterprises, pharmaceutical and health care, and scientific research institutes.



MODEL	DHG-9041A
Power	550W
Working Chamber Dimension L.xW.xH	350x350x350mm(43L)
Temperature Range	RT+10°C ~200°C
Temperature Resolution	0.1°C
Temperature Fluctuation	±1°C
Time	0-999min
Voltage	220v 50Hz