

## Automatic Upper Arm Blood Pressure Monitor

**M2 Basic**

**HEM-7120-E**



<b>Name</b>	M2 Basic
<b>Model</b>	HEM-7120-E
<b>Manufacturer</b>	OMRON HEALTHCARE Co., Ltd 53, Kunotsubo, Terado-cho, Muko, Kyoto 617-0002 Japan
<b>Manufacturing plants</b>	OMRON HEALTHCARE MANUFACTURING VIETNAM CO., LTD., Binh Duong Province, Vietnam
<b>Delivery set</b>	Device, cuff, storage case, blood pressure pass, batteries, instruction manual
<b>Measurement method</b>	Oscillometric
<b>Measurement range</b>	Pressure: 0mmHg to 299mmHg / Pulse: 40 to 180/min.
<b>Memory</b>	Last reading measurement
<b>Power</b>	4 “AAA” batteries 1.5V or AC/DC adapter (6V = 4W) (Adapter S-9515336-9, INPUT AC100-240V 50/60Hz 0.12A)