

Autoclave Sterilizer with Digital Display

100 L

LS-100HJ

Characteristics about steam sterilizer:

- 1.The steam sterilizer has fully stainless steel structure
- 2.Hand wheel type of quick-open door structure
- 3.Pressure safety lock
- 4.Indicator light indicates working state
- 5.Over temperature & over pressure auto-protection
- 6.Double scale indication pressure gauge
- 7.Self inflating type seal
- 8.Safe protection of water lacking
- 9.Automatically shut off with beep reminding after sterilization
- 10.Easy to operate safe and reliable
- 11.With two stainless steel sterilizing baskets



Technical Data:

Model	LS-35HJ	LS-50HJ	LS-75HJ	LS-100HJ
Technical data				
Chamber volume	35L(ϕ 318×450)mm	50L(ϕ 340×550)mm	75L(ϕ 400×600)mm	100L(ϕ 440×650)mm
working pressure	0.22MPa			0.14MPa
Working temperature	134°C			126°C
Max working pressure	0.23 Mpa			0.165MPa
Heat average	$\leq\pm 1^\circ\text{C}$			
Timer	0~60min			
Adjustment of temperature	105~134°C			105~126°C
Power	2.5Kw/AC220V.50Hz	3Kw /AC220V.50Hz	4.5KW /AC220V.50Hz	
Overall dimension	450×450×1010(mm)	510×470×1130(mm)	560×560×1120(mm)	540×560×1250(mm)
Transport dimension	570×550×1150(mm)	590×590×1280(mm)	650×630×1280(mm)	680×630×1370(mm)
G.W/N.W	72Kg/ 56Kg	88Kg/ 68Kg	100Kg/80Kg	110 Kg/85Kg