

DRY OVEN

DHG-9140A

- Microprocessor control.
- With double glass window.
- Digital display setting and actual temperature.
- Digital display running time.
- PID self-tuning, over temperature alarm.
- Automatic alarm after operation.
- Stainless-steel interior, chamber corners are semicircle.
- It is applicable to drying, baking, sterilization, dewaxing, solidifying for industrial and mining enterprises, pharmaceutical and health care, and scientific research institutes.



MODEL	DHG-9140A
Power	1500W
Working size(mm)	550x450x550
Temperature Range	50~300°C
Temperature Resolution	0.1°C
Temperature Fluctuation	±1°C
Time	0-999min
Voltage	110V 60Hz / 220v 50Hz
Dimensions(mm)	720x780x750
Weight	60kg