

Incubator RT-A19



Caracteristici generale

- I. Incubator uscat, compact.
- II. Pentru utilizarea cu analizorele semiautomate de biochimie
- III. Incubeaza tuburile la
 25°C / 30°C / 37°C / 45°C
- IV. Afisajul temperaturii cu LED-uri
- V. Control electronic cu microprocesor
- VI. Se poate utiliza cu 2 tipuri de tuburi, capacitatea maxima de incubare: 24 tuburi.

Caracteristici tehnice	
Modul de incalzire	Rezistor de metalina
Optiuni de temperatura	25°C / 30°C / 37°C / 45°C
Precizia temperaturii	$\pm 0.5^{\circ}\text{C}$
Direfenta de temperatura intre godeuri	$\pm 0.3^{\circ}\text{C}$
Viteza de incalzire	in domeniul $25\text{--}37^{\circ}\text{C}$ in aproximativ 6 minute (cu temperatura ambianta de 25°C)
Afisare temperatura	LED
Dimensiuni godeuri	20 cu $\Phi 13\text{mm}$ 、4 cu $\Phi 15\text{mm}$; adancime: 33mm
Putere element incalzire	70W
Dimensiuni L × I × H(mm)	272mm x 203mm x 188mm
Greutate neta	2.4kg
Alimentare	AC 85~265V 50~60HZ

in Romania prin:

BALMED s.r.l.

Bd. Burebista Nr. 1
Bl. D15, Sc.4, Ap. 125
031106 BUCURESTI
Tel.fax:021/327.52.69
021/327.09.21
office@balmed.ro
www.balmed.ro

Producator

Rayto Life and Analytical Sciences Co.,Ltd.
A/6F, Technology Building, Industry
Road, Shekou, Shenzhen P.R.China
Tel:+86 755 86168182 Fax:+86 755 86168796
E-mail: info@rayto.com Website: www.rayto.com

Rayto 雷杜