

PORTABLE AUTOCLAVE

AT-P24



Specification:

<i>PORTABLE AUTOCLAVE (DEEPENED)</i>	
<i>Art.No.</i>	<i>AT-P24</i>
<i>Sterilizing volume</i>	<i>24L</i>
<i>Working pressure</i>	<i>0.14-0.16Mpa</i>
<i>Working temperature</i>	<i>126°C</i>
<i>Max. safety pressure</i>	<i>0.165Mpa</i>
<i>Characteristics</i>	<i>Electric heated</i>