

# ***CENTRIFUGE***

## ***AVI-558***



### **Specification:**

<b><i>Maximum Rotor Speed</i></b>	<b><i>3500 RPM</i></b>
<b><i>Maximum Centrifuge Force</i></b>	<b><i>1200xg</i></b>
<b><i>Width</i></b>	<b><i>215 mm</i></b>
<b><i>Weight</i></b>	<b><i>4.7 kg ( Approx. Weight )</i></b>
<b><i>Connected Electrical Load</i></b>	<b><i>0.14 KVA</i></b>